Calbiochem®

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Laxminarayana R. Devireddy, Jose G. Teodoro, Fabien A. Richard, Michael R. Green(2001) Induction of Apoptosis by a Secreted Lipocalin That is Transcriptionally Regulated by IL-3 Deprivation. *Science* 293 (5531): 829-834

Nigel J. Waterhouse, Joshua C. Goldstein, Oliver von Ahsen, Martin Schuler, Donald D. Newmeyer, and Douglas R. Green(2001). Cytochrome c Maintains Mitochondrial Transmembrane Potential and ATP Generation after Outer Mitochondrial Membrane Permeabilization during the Apoptotic Process. *The Journal of Cell Biology*, 153(2): 319–328

Laxminarayana R. Devireddy, Jose G. Teodoro, Fabien A. Richard, Michael R. Green(2001). Induction of Apoptosis by a Secreted Lipocalin That is Transcriptionally Regulated by IL-3 Deprivation . *Science*, 293(5531): 829 – 834

Katja C. Zimmermann, Jean-Ehrland Ricci, Nathalie M. Droin, and Douglas R. Green(2002). The role of ARK in stress-induced apoptosis in Drosophila cells. *The Journal of Cell Biology*, 156, (6): 1077–1087

A Asumendi, MC Morales, A Alvarez, J Are´chaga and G Pe´rez-Yarza(2002). Implication of mitochondria-derived ROS and cardiolipin peroxidation in N-(4-hydroxyphenyl) retinamide-induced apoptosis. *British Journal of Cancer*, 86: 1951 – 1956

Caroline Sievers, Nick Platt, V. Hugh Perry, Michael P. Coleman, Laura Conforti (2003). Neurites undergoing Wallerian degeneration show an apoptotic-like process with annexin V positive staining and loss of mitochondrial membrane potential. *Neuroscience Research*, 46: 161-169

Vladimír Leksa, Samuel Godar, Herbert B. Schiller, Elke Fuertbauer, Arshad Muhammad, Katarina Slezakova, Vaclav Horejsi, Peter Steinlein⁵, Ulrich H. Weidle, Bernd R. Binder and Hannes Stockinger(2005). TGF- β -induced apoptosis in endothelial cells mediated by M6P/IGFII-R and mini-plasminogen. *Journal of Cell Science* , 118: 4577-4586

A Burchert, Y Wang, D Cai, N von Bubnoff, P Paschka, S Mu¨ller-Bru¨sselbach, OG Ottmann, J Duyster, A Hochhaus and A Neubauer (2005) . Compensatory PI3-kinase/Akt/mTor activation regulates imatinib resistance development. *Leukemia*, 19:1774–1782

Shawn D. Larson, Jing Li, Dai H. Chung, and B. Mark Evers(2007). Molecular mechanisms contributing to glutamine-mediated intestinal cell survival. *Am J Physiol Gastrointest Liver Physiol* 293: 1262–1271

ANNA-MARIA LOBASCIO, FRANCESCA-GIOIA KLINGER and MASSIMO DE FELICI (2007) . Isolation of apoptotic mouse fetal oocytes by AnnexinV assay. *Int. J. Dev. Biol.* 51: 157-160

Siriporn Jitkaew, Erika Witasp, Shouting Zhang, Valerian E. Kagan, and Bengt Fadeel(2009). Induction of caspase- and reactive oxygen speciesindependent phosphatidylserine externalization in primary human neutrophils: role in macrophage recognition and engulfment. *Journal of Leukocyte Biology*, 85: 427–437

Hyun Sun Cho, Seung Hee Chang, Youn Sun Chung, Ji Young Shin, Sung Jin Park, Eun Sun Lee, Soon Kyung Hwang, Jung Taek Kwon, Arash Minai Tehrani, Minah Woo, Mi Sook Noh, Huda Hanifah, Hua Jin, Cheng Xiong Xu, Myung Haing Cho(2009). Synergistic effect of ERK inhibition on tetrandrine-induced apoptosis in A549 human lung carcinoma cells. *J. Vet. Sci.* 10(1), 23-28

Christopher S Bryant, Sanjeev Kumar, Sreedhar Chamala, Jay Shah, Jagannath Pal, Mahdi Haider, Shelly Seward, Aamer M Qazi, Robert Morris, Assaad Semaan, Masood A Shammas, Christopher Steffes, Ravindra B Potti, Madhu Prasad, Donald W Weaver, Ramesh B Batchu(2010). Sulforaphane induces cell cycle arrest by protecting RB-E2F-1 complex in epithelial ovarian cancer cells. *Molecular Cancer*, 9:47

Jiamin Chen, Harriet E. Feilotter, Genevieve C. Pare´, Xiao Zhang, Joshua G.W. Pemberton, Cherif Garady, Dulcie Lai, Xiaolong Yang, and Victor A. Tron(2010). MicroRNA-193b Represses Cell Proliferation and Regulates Cyclin D1 in Melanoma. *The American Journal of Pathology*, 176(5):2520-2529