microRNA研究指南

AB Applied Biosystems

Ambion
THE RNA COMPANY®

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RecoverAll™ Total Nucleic Acid Isolation Kit

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
Figure 6. mirVana™ miRNA Isolation Kit for Efficient Recovery of miRNA. (A) Total RNA was isolated from the same mouse liver lysate using a double phenol/guanidinium extraction (DPGE) or the mirVana miRNA Isolation Kit procedure in triplicate (miRNA 1 to 3). The experiment was performed with two different mouse liver lysates. Each sample RNA (1 µg) was analyzed on a denaturing 15% polyacrylamide gel stained with ethidium bromide. (B) RNAs from the same gel were transferred to a membrane and probed for U2 snRNA and let-7 miRNA. The relative amount of small RNA in each lane was quantified with a phosphorimager. The graph shows the percentage of recovery with respect to the DPGE prep.

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
Figure 5. Differential Recovery of Small RNAs. (A) Total RNA was isolated from 1x10⁶ HeLa cells using three different techniques: monophasic phenol/chaotropic extraction (MPCE), binding on glass fiber filter in guanidinium solution (GFF) or double phenol/guanidinium extraction (DPGE). Purified RNA (1 μg) was resolved on a 1.2% denaturing agarose gel (top left panel) or 15% denaturing polyacrylamide gel (bottom left panel). The indicated mRNAs or small RNAs were detected by Northern blot or using the mirVana™ miRNA Detection Kit (right panels). (B) FirstChoice® Mouse Kidney Total RNA (20 μg) was either bound and eluted from a GFF or precipitated with 0.5 M NH₄OAc and 3 volumes of EtOH (EtOH PPT). The untreated or recovered RNAs (1 μg) were compared by gel analysis, Northern blot, or solution hybridization as described in (A).

Figure 7. Using the mirVana™ PARIST™ Kit to Analyze RNAi Effects. 1.5x10⁶ HUVEC cells were electroporated (800 V, 120 μs, 2 pulses, 0.5 s between pulses) with siRNA (10 μg) targeting GAPDH mRNA (GAPDH) or a negative control sequence (Neg) in 400 μl of siPORT™ siRNA Electroporation Buffer. Total RNA and protein were isolated with the mirVana PARIS Kit 48 hr after electroporation. One μg of total RNA was used to detect the indicated RNA species by Northern blot or with the mirVana miRNA Detection Kit. RNA probes were prepared by in vitro transcription (mRNA probes, MAXIsript® Kit) or by 5' end labeling (rRNA, miRNA, and siRNA probes, mirVana Probe & Marker Kit). Western blots were performed with 15 μg of total protein and antibodies specific for GAPDH or Ku p70 proteins.

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Figure 8. The flashPAGE™ Fractionator Apparatus.

“miRNA Certified” FirstChoice® Total RNA

Human	Total	RNA-Normal	CAT #	SIZE
Adrenal			AM7994	100 µg
Aorta			AM6844	100 µg
Left Atrium			AM6854	100 µg
Right Atrium			AM6858	100 µg
Bladder			AM7990	100 µg
Brain			AM7962	100 µg
Breast			AM7952	100 µg
Cervix			AM7992	100 µg
Colon			AM7986	100 µg
Distal Colon			AM6836	100 µg
Proximal Colon			AM6834	100 µg
Duodenum			AM6832	100 µg
Esophagus			AM6842	100 µg
Fallopian Tube			AM6862	50 µg
Heart			AM7966	100 µg
Ileum			AM6828	100 µg
Jujenum			AM6830	100 µg
Kidney			AM7976	100 µg
Liver			AM7960	100 µg
Lung			AM7968	100 µg
Lymph Node			AM7894	50 µg
Ovary			AM7974	100 µg
Pancreas			AM7954	100 µg
Pericardium			AM6852	100 µg
Placenta			AM7950	100 µg
Prostate			AM7988	100 µg
Skeletal Muscle			AM7982	100 µg
Small Intestine			AM7984	100 µg
Spleen			AM7970	100 µg
Stomach			AM7996	100 µg
Testes			AM7972	100 µg
Thymus			AM7964	100 µg
Thyroid			AM6872	50 µg
Trachea			AM6846	100 µg
Uterus			AM7892	100 µg
Vena Cava			AM6848	50 µg
Left Ventricle			AM6856	100 µg
Right Ventricle			AM6860	100 µg

ORDERING INFORMATION

Product	CAT #	SIZE
mirVana™ miRNA Isolation Kit	AM1560	40 purifications
mirVana™ PARIS Kit	AM1556	40 purifications
RecoverAll™ Total Nucleic Acid Isolation Kit for FFPE	AM1975	40 purifications


Figure 9. flashPAGE™ Fractionation of Decade Marker.


Ambion's Decade Marker, in a background of 10 µg mouse brain total RNA, was loaded onto a pre-cast flashPAGE Gel Cartridge and electrophoresed with the flashPAGE Fractionator. Two successive fractions were collected. For the first fraction, the lower buffer was collected and precipitated when the second dye indicator band had reached the lower end of the gel cartridge (middle lane). After adding fresh lower buffer to the apparatus, the sample was electrophoresed for 10 more minutes. The second lower buffer fraction was then collected and precipitated (right lane).

Mouse	Total	RNA-Normal	CAT #	SIZE
Liver			AM7810	200 µg
Brain			AM7812	200 µg
Thymus			AM7814	200 µg
Heart			AM7816	200 µg
Lung			AM7818	200 µg
Spleen			AM7820	200 µg
Testicle			AM7822	200 µg
Ovary			AM7824	200 µg
Kidney			AM7826	200 µg
Embryo (10–12 days)			AM7828	200 µg
Assorted			AM7800	10x25 µg

Rat Total RNA-Normal Tissue	CAT #	SIZE
Liver	AM7910	200 µg
Brain	AM7912	200 µg
Thymus	AM7914	200 µg
Heart	AM7916	200 µg
Lung	AM7918	200 µg
Spleen	AM7920	200 µg
Testicle	AM7922	200 µg
Ovary	AM7924	200 µg
Kidney	AM7926	200 µg
Embryo (10–12 days)	AM7928	200 µg
Assorted	AM7900	10x25 µg

Product	CAT #	SIZE
flashPAGE™ Fractionator	AM13100	1 apparatus
flashPAGE™ Pre-cast Gels (Type A)	AM10010	10 gels
flashPAGE™ Buffer Kit (Type A)	AM9015	1 kit
flashPAGE™ Reaction Clean-up Kit	AM12200	20 reactions

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

Performance Equivalent to Purified RNA

The TaqMan® MicroRNA Cells-to-CT™ Kit demonstrated equivalent performance vs. purified RNA when tested against a panel of 111 TaqMan® MicroRNA Assays. CT values from reactions using 1,000 HeLa cells prepared with TaqMan MicroRNA Cells-to-CT Kit (Y-axis) or purified total RNA (X-axis) as template were plotted. Results demonstrate similar expression patterns and equal sensitivity between both sample preparation procedures ($R^2 = 0.931$, slope = 0.977).

Linear Real-time RT-PCR using 10 to 100,000 Cells

A dilution series of 10–100,000 HeLa cells was processed in duplicate with the TaqMan® MicroRNA Cells-to-CT™ Kit, and with traditional purification. Two representative assays (Assay miR-16 #4373121, panel A, and Assay RNU48 #4373383, panel B) were performed in duplicate and averaged values were plotted as a function of cellular input.

ORDERING INFORMATION

Description	Size	Part Number
TaqMan® MicroRNA Cells-to-Ct™ Kit	100 lysis rxns	4391848
100 lysis reactions with gDNA removal	500 Real-Time PCR rxns	
200 miRNA cDNA synthesis reactions (15 µL)		
500 PCR (20 µL)		
TaqMan® MicroRNA Cells-to-Ct™ Kit	400 lysis rxns	4391996
400 lysis reactions with gDNA removal	2,000 Real-Time PCR rxns	
1,000 miRNA cDNA synthesis reactions (15 µL)		
2,000 PCR (20 µL)		

Please inquire about bulk pricing.

For Research Use Only. Not for use in diagnostic procedures.

RELATED PRODUCTS

Product Type	Description	Part Number
Accessories	MicroAmp® 96- & 384-Well Optical Adhesive Film	4311971
	RNaseZap® RNase Decontamination Solution	AM9782
	Nuclease-free Water (not DEPC-treated)	AM9938
	RT-PCR Grade Water	AM9935
Assays	TaqMan® MicroRNA Assays	
	A comprehensive collection of >900 pre-designed, pre-optimized probe and primer sets for quantitative gene expression analysis of human, mouse, rat, <i>Arabidopsis</i> , <i>Drosophila</i> , and <i>C. elegans</i> mature miRNAs http://mirna.appliedbiosystems.com	
Instruments	StepOne™ Real-Time PCR System	4376357
	StepOnePlus™ Real-Time PCR System	4376600
	Applied Biosystems 7300 Real-Time PCR System	4351101
	Applied Biosystems 7500 Real-Time PCR System	4351104
	Applied Biosystems 7900HT Fast Real-Time PCR System, Standard 96-Well Block	4329003
	Applied Biosystems 7900HT Fast Real-Time PCR System, Standard 384-Well Block	4329001
RT	TaqMan® MicroRNA Reverse Transcription Kit	200 rxns 4366596
Master Mixes	TaqMan® Universal PCR Master Mix	200 rxns 4324018
		400 rxns 4364341
		1,000 rxns 4364343
		2,000 rxns 4326614
Plates	MicroAmp® Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL	4366932
	MicroAmp® Optical 384-Well Reaction Plate	4343370
	MicroAmp® Optical 96-Well Reaction Plate	4316813

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TaqMan MicroRNA Assays——实时定量PCR的革命

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
④ Dilute cDNA & Load each TaqMan® Array

Diluted RT product and TaqMan® Universal PCR Master Mix

10 min

⑤ Run TaqMan Array


Applied Biosystems 7900HT Fast Real-Time PCR System

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TaqMan MicroRNA Assay Endogenous Controls

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TaqMan MicroRNA Reverse Transcription (RT) Kit

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ORDERING INFORMATION

Content		Input Quantity	Species	Megaplex™ Pools P/N	TaqMan™ Array P/N	TaqMan™ PreAmp Master Mix P/N	TaqMan™ MicroRNA RT Kit P/N	TaqMan™ Universal PCR Master Mix, No AmpErase® UNG P/N
A pool only	Unlimited Sample (without PreAmp Primers)	Human	4399966	4398965	N/A	4366596 or 4366597	4324018	
		Rodent	4399970	4398967				
		Limited Sample (includes PreAmp Primers)	Human	4401009	4398965			4391128
			Rodent	4401090	4398967			
	B pool only	Unlimited Sample (without PreAmp Primers)	Human	4399968	4398966			N/A
			Rodent	4399972	4398968			
		Limited Sample (includes PreAmp Primers)	Human	4401010	4398966			4391128
			Rodent	4401011	4398968			
A and B pools	Unlimited Sample (without PreAmp Primers)	Human	4400928	4400238	N/A			
		Rodent	4400926	4400239				
	Limited Sample (includes PreAmp Primers)	Human	4401091	4400238	4391128			
		Rodent	4401012	4400239				

Additional TaqMan MicroRNA Assay Products

	P/N	Contents
TaqMan™ MicroRNA Assays (Individual)	Multiple ‡	<ul style="list-style-type: none"> • RT primer tube: 5X reaction mix • TaqMan™ Assay tube: 20X reaction mix • 50 RT (15µL) and 150 real-time PCR reactions (20µL)
TaqMan™ MicroRNA Endogenous Controls	Multiple ‡	<ul style="list-style-type: none"> • RT primer tube: 5X reaction mix • TaqMan Assay tube: 20X reaction mix • 50 RT (15µL) and 150 real-time PCR reactions (20µL)
TaqMan™ MicroRNA Assay Sets	Multiple ‡	Per Assay: <ul style="list-style-type: none"> • RT primer tube: 5X reaction mix • TaqMan Assay tube: 20X reaction mix • 50 RT (15µL) and 150 real-time PCR reactions (20µL)
TaqMan™ MicroRNA RT Kit		<ul style="list-style-type: none"> • 10X RT buffer • dNTP mix w/dTTP (100mM total) • RNase Inhibitor (20U/µL) • MultiScribe™ RT enzyme (50U/µL)
200 reactions (100 rxns without PreAmp, 200 with PreAmp)	4366596	
1,000 reactions (500 rxns without PreAmp, 1,000 with PreAmp)	4366597	

‡ For the latest information on TaqMan MicroRNA Assays, including a list of available assays, visit mirna.appliedbiosystems.com

Related Products Required for Use with TaqMan MicroRNA Assays

TaqMan™ Universal PCR Master Mix, No AmpErase™ UNG

TaqMan™ PreAmp Master Mix

mirVana™ miRNA Isolation Kit

Other Related Products for Use with TaqMan MicroRNA Assays

TaqMan™ MicroRNA Cells-to-CT™ Kit

Anti-miR™ miRNA Molecules, Controls, and Libraries

Pre-miR™ miRNA Molecules, Controls, and Libraries

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
Figure 22. Activity of Pre-miR™ miRNA Precursor Molecules and Anti-miR™ miRNA Inhibitors.

ORDERING INFORMATION		
Product	CAT #	SIZE
mirVana™ miRNA Probe and Marker Kit	AM1554	30 rxns; 10 marker rxns
Anti-miR™ miRNA Inhibitors	AM17000	5 nmol
	AM17001	20 nmol (4x5 nmol)
order online at www.ambion.com/catalog/mirna_search.php		
Anti-miR™ miRNA Inhibitors – custom defined sequence	AM17003	20 nmol
order online at www.ambion.com/catalog/mirna_search.php		
Anti-miR™ Negative Control #1	AM17010	5 nmol

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ORDERING INFORMATION

Product	CAT #	SIZE
Pre-miR™ miRNA Precursor Molecules	AM17100	5 nmol
order online at www.ambion.com/catalog/mirna_search.php	AM17101	20 nmol (4x5 nmol)
Pre-miR™ miRNA Precursor Molecules custom defined sequence	AM17103	20 nmol
order online at www.ambion.com/catalog/mirna_search.php		
Pre-miR™ Negative Control #1	AM17110	5 nmol
Pre-miR™ Negative Control #2	AM17111	5 nmol

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ORDERING INFORMATION

Product	CAT #	SIZE
siPORT™ NeoFX™ Transfection Agent	AM4510	0.4ml
	AM4511	1.0ml

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