

Easy assay setup and high-quality data using the LiquiChip™ System

The LiquiChip™ Protein Suspension Array System is a bead-based platform that enables a wide range of protein interaction arrays to be carried out quickly and with high sensitivity (see *QIAGEN News* 2002 No 2, 1). Here we show how assays using 6xHis-tagged proteins can be quickly and easily developed and used to obtain high-quality assay data.

Easy assay setup

One of the major advantages of the LiquiChip System is that it enables QIAexpress® 6xHis technology to be used in xMAP™ bead-based assays. This means that a single chemistry can be used for all steps of the protein-handling process, from purification to immobilization and assay. Using a common immobilization protocol for 6xHis-tagged proteins increases standardization in assays and eliminates the need for lengthy and cost-intensive optimization trials (Figure 1). Proteins are bound to beads in a batch procedure and in most cases can be added directly to LiquiChip Assays. Table 1 shows a comparison of the steps involved in coupling 6xHis-tagged proteins to LiquiChip Ni-NTA and Penta-His™ Beads with those required for covalent protein immobilization.

In addition to being quicker, immobilization of 6xHis-tagged proteins on Ni-NTA and Penta-His Beads requires less protein per assay point. Proteins are bound in a directed manner, maximizing accessibility to interaction partners and increasing signal intensity.

High-quality data sets

Assay data quality is reflected in signal-to-noise and signal-to-background ratios. The Z-factor (see Table 2 and reference 1) is a dimensionless, simple statistical characteristic and the most suitable parameter for assay quality assessment. The Z-factor reflects assay signal dynamic range and variation associated with signal measurements.

To determine the Z-factor of a typical LiquiChip assay, a hypothetical drug-screening assay using proteases was established (Figure 2). 6xHis-tagged Thioredoxin-Tag-100 was immobilized on LiquiChip Ni-NTA Beads, incubated with differing amounts of enterokinase, and residual noncleaved protein was quantified after 3 hours incubation using Tag-100 specific antibodies. Data were collected from three processed 96-well plates ($n = 288$). ▶

Reference

1. Zhang, J.H., Chung, T.D., and Oldenburg, K.R. (1999) A simple statistical parameter for use in evaluation and validation of high-throughput screening assays. *J. Biomol. Screen.* **4**, 67.

Table 1. Steps involved in coupling 6xHis-tagged proteins to LiquiChip Ni-NTA and Penta-His Beads in comparison with covalent protein immobilization

	Ni-NTA Bead coupling procedure	Penta-His Bead coupling procedure	Covalent protein coupling procedure
Volume of beads required for 500 assay points	250 µl	250 µl	1 ml
Washing step	Not required	Not required	2 centrifugation steps
Activation step	Not required	Not required	EDC/Sulfo-NHS (1 h)
Protein required for 500 assay points	<5 µg	<5 µg	25–250 µg (optimization recommended)
Incubation	1 h or overnight	1 h or overnight	0.5 h – overnight
Additional steps required	Ready to use	Ready to use	Further wash step Enumeration of beads

Table 2. The Z-factor

The Z-factor is based on assay signal-to-noise and signal-to-background ratios.

$$Z\text{-factor} = 1 - \frac{3\sigma \text{ positive control} + 3\sigma \text{ negative control}}{\text{mean positive control} - \text{mean negative control}}$$

For the experiment presented in Figure 2 the following values were calculated.

Background (negative control): Mean = 18, σ (standard deviation) = 2.5
 Positive control: Mean = 2524, σ (standard deviation) = 66.5

These data give a Z-factor for the assay shown in Figure 2 of 0.92.

Sensitive Multiplex Assays




Figure 1 Multiplex ELISA assay in which three different proteins were bound to LiquiChip Ni-NTA Beads using a common procedure, and used to quantify protein-specific antibodies.

Highly Reproducible Assays




Figure 2 High signal-to-background ratio and low background noise result in high Z-values. As a positive control, protease activity was inhibited with the serine-protease inhibitor aprotinin (10 ng/well). Data derived from 3 microtiter plates, processed on three successive days, are shown.

By enabling fast assay development and setup and delivering high-quality data, the LiquiChip System is a valuable addition to any lab performing protein interaction assays. For more information on the

LiquiChip System and how it can help you to streamline assay development visit www.qiagen.com or call QIAGEN Technical Services. ■