High Pure PCR Cleanup Micro Kit 创新灵活的核酸纯化产品

Nina Lassonczyk, Constanze Stadler, Jenna Liedtke 和 Horst Donner*

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Table 1: Application selection guide.				
For 100 µl liquid sample or 100 mg agarose gel		Option 1	Option 2	Option 3
		200 µl binding buffer+200 µl binding enhancer(40%)	300 µl binding buffer+100 µl binding enhancer(20%)	400 µl binding buffer
Purification of liquid sample,100 µl	Laqbeling or other reaction products 100 bp to 5 kb	+		
	PCR products 100 bp to 5 kb	+		
	DNA fragments for sequencing			
Purification from agarose gel,100 mg	DNA Fragments 100 bp to 5 kb		+	+
Removal of low molecular weight DNA	Primer up to 25 bases			
	Primer-dimer up to 70 bp		+	+

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

Figure 2: Agarose gel electrophoresis and ethidiumbromide staining of 341 bp PCR products.

Table 2: Yield of PCR fragments isolated with increasing amounts of binding enhancer.

Binding conditions	260/280 ratio	Yield (ug)
No binding enhancer (option 3)	1.9	0.6
10% binding enhancer	1.9	0.9
20% binding enhancer	1.8	1.2
40% binding enhancer	1.8	1.2


Figure 3: Melting curve analysis after the last PCR cycle using a LightCycler 1.5 Instrument.

Table 3: Purification of an 1,100 bp cDNA fragment in combination with depletion of Cy5 dUTP.

Sample ID	Yield (ug)	260/280 ratio	Cy5 dye (mmol)
Supplier Q	0.1	1.9	not detectable
Supplier M	1.1	2	0.22
High Pure PCR Cleanup Micro Kit + 20% binding enhancer (option 2)	1.8	2.1	not detectable
High Pure PCR Cleanup Micro Kit + 40% binding enhancer (option 1)	1.9	2.1	0.01

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www.ras-biosino.com

FuGENE® HD 转染试剂

**Measure the results of your transfection,
not your transfection reagent.**


Figure 1. Minimize off-target effects by using FuGENE® HD Transfection Reagent.
FuGENE® HD Transfection Reagent or a reagent from another supplier (L2K) was used to transfect MCF-7 cells (RXCC® HTB-22®). Subsequent microarray expression profiling experiments demonstrated that L2K significantly altered the expression levels of six times more genes than FuGENE® HD Transfection Reagent. (View the complete article online in *Biochemistry* (2006) 4 at www.roche-applied-science.com/publications/Biochemistry.htm)

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