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
■ PrimeSTAR GXL DNA Polymerase相关参考文献

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癌症 | Illumina测序

Sakaguchi, H., *et al.* Exome sequencing identifies secondary mutations of *SETBP1* and *JAK3* in juvenile myelomonocytic leukemia. *Nature Genetics* **45**, 937–941 (2013).

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Kakizaki, E., *et al.* Detection of diverse aquatic microbes in blood and organs of drowning victims: First metagenomic approach using high-throughput 454-pyrosequencing. *Forensic Sci.* **220**, 135–146 (2012).

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Hosomichi, K., *et al.* Phase-defined complete sequencing of the HLA genes by Next-Generation Sequencing. *BMC Genomics* **14**, 355 (2013).

Izawa, K., *et al.* Detection of base substitution-type somatic mosaicism of the NLRP3 gene with >99.9% statistical confidence by massively parallel sequencing. *DNA Res.* **19**, 143–152 (2012).

Shiina, T., *et al.* Super high resolution for single molecule-sequence-based typing of classical HLA loci at the 8-digit level using next generation sequencers. *Tissue Antigens* **80**, 305–316 (2012).

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Al-Khedery, B., *et al.* Structure of the type IV secretion system in different strains of *Anaplasma phagocytophilum*. *BMC Genomics* **13**, 678 (2012).

Novitsky, V., *et al.* Long-range HIV genotyping using viral RNA and proviral DNA for analysis of HIV drug resistance and HIV clustering. *J. Clin. Microbiol.* **53**, 2581–2592 (2015).

Ode, H., *et al.* Quasispecies analyses of the HIV-1 Near-full-length genome with illumina MiSeq. *Front. Microbiol.* **6**, 1258 (2015).

Takeda, T. S., Sugiyama, M. A. & Kanto, T. A. Establishment of a practical workflow using a next generation sequencer to search for novel SNPs associated with the anti – HCV treatment response. *Juntendo Med. J.* **60**, 568–575 (2014).

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Sugiura, M., *et al.* Induced pluripotent stem cell generation-associated point mutations arise during the initial stages of the conversion of these cells. *Stem Cell Reports* **2**, 52–63 (2013).

*also used Titanium Taq

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Hosomichi K. *et al.* Phase-defined complete sequencing of the HLA genes by next-generation sequencing. *BMC Genomics* **28**, 355 (2013).

Hosomichi, K., Shiina, T., Tajima, A. & Inoue, I. The impact of next-generation sequencing technologies on HLA research. *J. Hum. Genet.* **60**, 665–73 (2015).

Mayor, N. P. *et al.* HLA typing for the next generation. *PLoS One* **10**, (2015). Shiina, T. *et al.* Rapid evolution of major histocompatibility complex class I genes in primates generates new disease alleles in humans via hitchhiking diversity. *Genetics* **173**, 1555–1570 (2006).

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Ehrenberg, P. K. *et al.* High-throughput multiplex HLA genotyping by next-generation sequencing using multilocus individual tagging. *BMC Genomics* **15**, 864 (2014).

Ozaki, Y., *et al.* HLA-DRB1, -DRB3, -DRB4 and -DRB5 genotyping at a super-high resolution level by long range PCR and high-throughput sequencing. *Tissue Antigens* **83**, 10–16 (2014).

Superior amplification of GC-rich templates with PrimeSTAR GXL DNA Polymerase

Data kindly provided by: Alexander Ebralidze, Scientist, Beth Israel Deaconess Medical Center

Introduction

PCR amplification of GC-rich templates is often hindered by formation of secondary structures and the requirement for high melting temperatures. **PrimeSTAR GXL DNA Polymerase** is a high-fidelity PCR enzyme that allows efficient amplification of the most challenging templates. With this enzyme, amplification of GC-rich templates is possible without reaction optimization or inclusion of additives such as DMSO or Betaine.

In this experiment, PrimeSTAR GXL DNA Polymerase was compared to other commercially available enzymes for the amplification of a target with high GC content. Only the PrimeSTAR GXL enzyme was able to efficiently amplify the GC-rich template without additives or special reaction conditions.

Results

Under standard reaction conditions, only PrimeSTAR GXL DNA Polymerase produced a single amplification product of the expected size (Figure 1). PCR product was successfully obtained with the competitor enzymes only when the reactions included DMSO (5% final) and Betaine (1 M final) (Figure 2, left). Very little product was obtained with PrimeSTAR GXL polymerase in the presence of these additives (Figure 2, right).


Figure 1. Amplification products obtained under standard reaction conditions. The expected product size is 603 bp. Pfu: Pfu polymerase; Co. A: Company A master mix; Co. D: Company D master mix; PS: PrimeSTAR GXL polymerase.


Figure 2. Amplification products obtained from reactions containing DMSO and Betaine. The expected product size is 603 bp. Pfu: Pfu polymerase; Co. A: Company A master mix; Co. D: Company D master mix; PS: PrimeSTAR GXL polymerase.

Conclusions

In this experiment, PrimeSTAR GXL DNA Polymerase was the only PCR enzyme able to amplify a sequence from template DNA with a GC content >75% without additives. PrimeSTAR GXL DNA polymerase was able to amplify this target without special reaction conditions (a standard two-step PCR protocol was used) and without additives. In contrast, all of the other competitor enzymes tested were only able to amplify when DMSO and Betaine were added to the reaction. These results indicate that **PrimeSTAR GXL DNA Polymerase** is ideal for optimization-free PCR amplification of GC-rich targets.

Methods

A 603-bp region of the *CEBPA* gene promoter that contains 75.29% GC content was amplified from human genomic DNA (isolated from the K562 cell line) using either PrimeSTAR GXL DNA Polymerase, a master mix from Company A, a master mix from Company D, or *Pfu* DNA polymerase. Reactions were assembled according to the manufacturer's recommendations. All reactions were run on a MultiGene Thermal Cycler (Labnet) following the cycling conditions in Table 1 below. Amplified products were analyzed by gel electrophoresis.

Table 1. PCR cycling conditions					
PrimeSTAR GXL polymerase (2-step PCR)			Company A, Company D, and <i>Pfu</i> polymerase (3-step PCR)		
98°C	10 sec	35 cycles	95°C	30 sec	35 cycles
68°C	20 sec		60°C	15 sec	
			72°C	30 sec	

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
■ PrimeSTAR Max DNA Polymerase相关参考文献

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Tanii, H., *et al.* Induction of Cytochrome P450 2A6 by bilirubin in human hepatocytes. *Pharmacol. Pharmacy* **4**, 182–190 (2013).

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Yamamoto, K., *et al.* A novel TRB@/NOTCH1 fusion gene in T-cell lymphoblastic lymphoma with t(7;9) (q34;q34). *European J. Haematology* **90**, 68–75 (2012).

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Yamanaka, K., *et al.* Rapid detection for primary screening of influenza A virus: microfluidic RT-PCR chip and electrochemical DNA sensor. *Analyst* **136**, 2064–2068 (2011).

■ PrimeSTAR Max DNA polymerase相关Tech note

Go from DNA to data with a fast polymerase that's also astonishingly accurate


PrimeSTAR Max DNA Polymerase

- Highest fidelity of any commercially available polymerase
- Fastest extension speed available
- Formulated as a hot-start premix with dNTPs and buffer: just add primers, DNA, bring to final volume with dH₂O, and go.

Introduction

A fast polymerase for high fidelity PCR—save your sanity as well as your time

Let's face it, there's nothing more frustrating than a failed PCR. When all that stands between you and understanding your gene of interest is DNA amplification, every moment counts. Rely on **PrimeSTAR Max DNA Polymerase** for cloning and expression, structure-function, and genetic variation studies.

Using a fast polymerase to detect influenza. Scientists at Osaka University developed a microfluidic device for rapid detection of influenza A (shown in color in the figure above). Using a mix of PrimeScript RTase, PrimeSTAR Max polymerase, and SpeedSTAR polymerase, RT-PCR was completed in 15 minutes.

Citation: Yamanaka, K. *et al.* Rapid detection for primary screening of influenza A virus: microfluidic RT-PCR chip and electrochemical DNA sensor. *Analyst* **136**, 2064–8 (2011).

Results

Experimental example: PrimeSTAR Max vs. other commercially available high-fidelity polymerases using a 10-second extension step

Human genomic DNA was used as template for fast PCR. Products 0.5 to 7.5 kb long were amplified (extension time, 10 sec) using PrimeSTAR Max enzyme or polymerases from three other manufacturers. High yield was observed when PrimeSTAR Max polymerase was used for products up to 4 kb in length.


Figure 1. Fast PCR amplification of human genomic DNA using PrimeSTAR Max DNA Polymerase or high-fidelity polymerases provided by other manufacturers.

Compare extension-step recommendations

Is every high-fidelity polymerase capable of fast extension? No. Many high-fidelity polymerases on the market lack sufficient extension speed and are incompatible with fast cycling conditions. Refer to the table below for recommended extension step times for PrimeSTAR Max DNA Polymerase and three other commercially available high-fidelity enzymes.

Polymerase	Total extension step time per cycle for a 4-kb product
PrimeSTAR Max polymerase	10 sec minimum*
Company T polymerase	120 sec
Company N polymerase	120 sec
Company I polymerase	240 sec

*See data in Figure 1 for an example of amplification of a 4-kb product with a 10 sec extension step. Routine recommendations: when template quantity is 200 ng/50 μ l or less, an extension step of 5 sec/kb is recommended. When template quantity exceeds 200 ng/50 μ l, an extension step of 30–60 sec/kb is recommended.

Conclusions

PrimeSTAR Max DNA Polymerase offers the highest fidelity and fastest extension speed of any commercially available polymerase. It was shown to provide superior yields for a series of amplification products up to 4 kb in length when compared to polymerases from other manufacturers.

PrimeSTAR Max DNA Polymerase also requires a far shorter extension time than other commercially available high-fidelity polymerases, making it a better choice for high-fidelity PCR performed under fast cycling conditions.

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